

# PhotoGenica SV

## Pulsed Dye Surgical Laser

Changing the management of laryngeal disease by targeting the lesion's vasculature



PhotoGenica<sup>®</sup> SV



- PDL non-ablative vascular specificity provides selective coagulation of micro-vasculature
- Improved hemostasis
- Non-contact treatment minimizes trauma to surrounding tissue
- Flexible fiber optic delivery for superior access and control



10 Elizabeth Drive, Chelmsford, MA 01824 USA • Tel. (978) 256-4200  
Toll Free: (800) 886-2966 Fax (978) 256-6556  
[www.cynosurelaser.com](http://www.cynosurelaser.com)

# Specifications

Wavelength (nm) .....	585
Pulse Width (msec) .....	0.45
Treatment Energy Range (J) .....	0.4 - 2.2
Rep. Rate (Hz) .....	2
Aiming Beam .....	635 nm diode laser
Electrical Requirements .....	120V/20A or 220V/20A
Size (HxWxD) .....	15" x 22" x 33" (38cm x 55cm x 84cm)
Weight (lbs).....	175
Dye Kit (delivered pulses) .....	20,000

PhotoGenica<sup>®</sup> SV  
Surgical Dye Laser



10 Elizabeth Drive, Chelmsford, MA 01824 USA • Tel. (978) 256-4200  
Toll Free: (800) 886-2966 Fax (978) 256-6556  
[www.cynosurelaser.com](http://www.cynosurelaser.com)